



คุณสมบัติของเครื่องผลิตออกซิเจน ยี่ห้อ OWGELS รุ่น OZ-301TWO

Model Parameter:

Project	Parameter
Weight	15kg
External dimension	300X260X520mm
Power	250W
Average sound pressure level	≤50dB (A)
Oxygen pressure	30-60KPa
Oxygen flow	3L/min±10%
Oxygen concentration	≥90%

Model Characteristics:

Medical oxygen supply: oxygen concentration is up to 93%.

Efficient and stable: produce oxygen for 24 hours.

One machine with multi-functions: atomization and oxygen concentration detection, meet more demands