



คุณสมบัติของเครื่องผลิตออกซิเจนขนาด 5 ลิตร ยี่ห้อ OWGELS รุ่น OZ-501TWO

Model Parameter:

Project	Parameter
Weight	19KG
External dimension	340X300X650mm
Power	350W
Average sound pressure level	$\leq 50\text{dB}$ (A)
Oxygen pressure	30-60kpa
Oxygen flow	5L/ min $\pm 10\%$
Oxygen concentration	$\geq 90\%$

Model Characteristics:

Imported originals: German oil-free compressor, French CECA molecular sieve to ensure high quality.

Efficient oxygen production: flow rate of 5L/min and high concentration of 93% $\pm 3\%$ to really meet medical demand.

Stable oxygen supply: continue to produce oxygen for 24 hours and operate for 365 days without stopping to guard healthy at all time. Intelligent alarm: AB-level intelligent alarm system to let it safer. Atomization function: fine atomization and effective treatment. Patent protection: 8 EU-certified patented technologies, authoritative and reliable.