

ห้างหุ้นส่วนจำกัด ซิกม่า ไบโอเทค

6/7 หมู่บ้านสุขสันต์ 9 ซ.เพชรเกษม 110 แขวงหนองค้างพลู เขตหนองแขม กทม.10160

โทร. 02-147-4608 , 094-549-8211

WWW.ZIGMABIOTECH.COM , Line ID : ZIGMABIOTECH



เครื่อง Hydrocolator ยี่ห้อ XIANGYU รุ่น XY-SRF-II

คุณสมบัติ

1. Hot wet compress device's rated voltage: 220V.
2. Hot wet compress device's rated operating frequency: 50Hz.
3. Hot wet compress device's input normal rated power: 1500W.
4. Volume of water tank: 80L.
5. Weight of hot wet damp compress device with water and hot compress bag: 150Kg.
6. Temperature range: adjustable within 0~99°C.
7. Stop heating over temperature range: 98°C-100°C.
8. Heating time: it needs 180 min to heat up from room temperature 25°C to 70°C.

9. Cooling time: it needs 180 min to cool down from 70°C to 25°C.

10. Related size of the hot wet compress device

10. 1 External dimensions of constant temperature heating box: 53cm (W) × 81cm (L) × 91cm (H) ±3cm

10. 2 Dimensions of inner box of the constant temperature heating box: 40cm (W) × 41cm (L) × 49cm (H)

±3cm.

11. Normal working conditions

12. 1 Environment temperature range: -18°C~60°C.

13. 2 Relative humidity range: 10%~100%.

14. 3 Atmospheric pressure: 50kpa~160kpa

15. 4 Power supply: a.c220V±10%; 50Hz±2%