

ห้างหุ้นส่วนจำกัด ซิกม่า ไบโอเทค

6/7 หมู่บ้านสุขสันต์ 9 ซ.เพชรเกษม 110 แขวงหนองค้างพลู เขตหนองแขม กทม.10160

โทร. 02-147-4608 , 094-549-8211

WWW.ZIGMABIOTECH.COM , [Line ID : ZIGMABIOTECH](#)



เครื่อง Hydrocolator ยี่ห้อ XIANGYU รุ่น XY-SRF-I

คุณสมบัติ

1. 13 bags in total:10 standard bags, 1 longer bag, 1 bag for neck, 1 bag for hand,
2. Intelligence operation:device will start to work /stop to work according to your previous setting
- 3, Duplicate protection against over heated and alarm will sound as the water lower than the required water volume
- 4, Constant temperature control: the standard indoor temperature to 99° adjustable
5. 220V±10% 50Hz±1Hz
6. Size:720*470*1030mm
7. Weight: device +all bags:85kg

8, Volume of water tank: 70L

9. Heating time: it needs 180 min to heat up from 25°C to 70°C

10. Cooling time: it needs 180 min to cool down from 70°C to 25°C

11. Size of Hot compress bags :Hand hot packs: 17cm×32cm±3cm, Standard hot packs:

25cm×30cm±3cm, Neck and shoulder hot packs: the length is

61cm±3cm, Bigger hot packs: 25cm×61cm±3cm.