



คุณสมบัติของเครื่องผลิตออกซิเจนขนาด 1 ลิตร(ขนาดพกพา)
ยี่ห้อ OWGELS รุ่น ZY-304

Model Parameter:

Project	Parameter
Weight	6.1kg
External dimension	350X235X280mm
Power	125W
Average sound pressure level	$\leq 50\text{dB}$ (A)
Oxygen pressure	30-60KPa
Oxygen flow	1L/min $\pm$ 10%
Oxygen concentration	$\geq 90\%$

Model Characteristics:

With an oxygen concentration of up to 93%, the oxygen is pure Child lock Ultra-silence, HD screen, portable and easy to carry